

Product Datasheet

Ephrin B1 (phospho-Tyr317) antibody (orb183417)

Description	Ephrin B1 (phospho-Tyr317) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Ephrin B1 around the phosphorylation site of Tyr317 (PH(p-Y)EK)
Target	phospho-Ephrin B1 (Tyr317)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34 kDa
Uniprot ID	P98172
Expiration Date	12 months from date of receipt.