

Product Datasheet

CARD16 antibody (orb18343)

Description

Goat polyclonal antibody to CARD16



Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, IF

Applications

Target

COP1 / PSEUDO-ICE (Isoform 2)

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage


Aliquot and store at -20°C. Minimize freezing and thawing.

Note


For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:4000. WB: Preliminary experiments gave an approx 28-30kDa band in human colon, pancreas, placenta and tonsil lysates after 0.3 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 10.7kDa according to NP_443121.1. The 28-30kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported



Immunofluorescence analysis of paraforma...



Flow cytometric analysis of paraformalde...

Protein Sequence

ETLGLSAGPIPGN

Clonality

Polyclonal

MW

10.7

Entrez

114769

Dilution Range

ELISA: 1:4000, IHC-P: 10ug/ml

Expiration Date

12 months from date of receipt.