


Product Datasheet

MTMR9 antibody (orb18344)

Description	Goat polyclonal antibody to MTMR9
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC
Target	MTMR9 (C Terminus)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.IHC: Immunohistochemistry: This product has been successfully used on Rat Hypothalamus (Yanagiya et al, Hum Mol Genet. 2007 Dec15;16(24):3017-26. PMID: 17855449).WB: No signal obtained yet but low background observed in Human Brain and Kidney lysates at up to 1 µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	AELETEDGMQESP
Clonality	Polyclonal
MW	63.5
Entrez	66036
Dilution Range	ELISA: 1:8000, WB: 1 µg/ml, IHC-P: 5-10µg/ml
Expiration Date	12 months from date of receipt.



orb18344
(5ug/ml)
staining of
paraffin e...