

Product Datasheet

ELMO2 antibody (orb18347)

Description	Goat polyclonal antibody to ELMO2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Target	ELMO2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 80kDa band seen in Human Brain lysate. Recommended for use at 1-3 µg/ml.
Protein Sequence	PKEPSSYDFVYHYG
Clonality	Polyclonal
MW	82.6; 72.8
Entrez	63916
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 3-6 µg/ml
Expiration Date	12 months from date of receipt.



orb18347
staining
(2ug/ml) of
Human Brai...