


Product Datasheet

CBL antibody (orb18349)

Description	Goat polyclonal antibody to CBL
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CBL
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 120kDa band seen in 3T3 lysate. Recommended for use at 1-3 µg/ml.
Protein Sequence	REFVSISSPAHVAT
Clonality	Polyclonal
MW	104
Entrez	867
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb18349
staining
(2ug/ml) of
3T3 lysate...