



Product Datasheet

COMMD8 antibody (orb1835)


Description	COMMD8 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human COMMD8 (101-183/183 aa)
Target	COMMD8
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21 kDa
Uniprot ID	Q9NX08
Expiration Date	12 months from date of receipt.



Western blot analysis of Eye (Mouse) Lys...



Immunohistochemical staining of COMMD8 a...



Immunohistochemical staining of COMMD8 a...