

Product Datasheet

RACGAP1 antibody (orb18360)

Description

Goat polyclonal antibody to RACGAP1

Species/Host

Goat

Reactivity

Bovine, Canine, Human

Conjugation

Unconjugated

Tested

ELISA, FC, IF, IHC, WB

Applications

Target

RACGAP1 / MgcRacGAP

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded Human Testis, shows nuclear staining of many cells. Recommended concentration, 1µg/ml.WB: Approx 60kDa band observed in Human Testis lysates (calculated MW of 71.0kDa according to NP_037409.2). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Toure et al., J Biol Chem 1998 Mar 13;273(11):6019-23; PMID: 9497316). Recommended concentration: 0.1-0.3 µg/ml.

Protein Sequence

GRQGNFFASPMKL

Clonality

Polyclonal

MW

71.0; 63.9; 72.3; 61.9

Entrez

29127

Dilution Range


ELISA: 1:32000, WB: 0.3-1 µg/ml, IHC-P: 1 µg/ml

Expiration Date


12 months from date of receipt.



Western blot analysis of Human Testis ly...



Immunohistochemical staining of Human Te...



Western blot analysis of HEK293 lysate u...