



Product Datasheet

FAM38A antibody (orb183601)


Description	FAM38A antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human FAM38A (2451-2521/2521 aa)
Target	PIEZO1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	287 kDa
Uniprot ID	Q92508
Dilution Range	IHC-P=1:100-500, ICC=1:100, IF=1:100-500, Flow-Cyt=1ug/Test, IHC-F=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Paraformaldehyde-fixed, paraffin embedde...



Paraformaldehyde-fixed, paraffin embedde...



IF analysis of HeLa cell lysate using FA...