



Product Datasheet

ELMO3 antibody (orb18361)

Description	Goat polyclonal to ELMO3
Species/Host	Goat
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ELMO3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 70kDa band seen in Human Brain lysate. Recommended for use at 1-3 µg/ml.
Protein Sequence	TNFNFCYDCSIAEP
Clonality	Polyclonal
MW	70
Entrez	79767
Dilution Range	ELISA: 1:16000, WB: 0.1-0.3 µg/ml, IHC-P: 5-10µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Brain lys...



Immunohistochemical staining of Human Pr...