


Product Datasheet

HAT4 antibody (orb183707)

Description	HAT4 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human HAT4/FLJ14154 (1-100/242 aa)
Target	HAT4/FLJ14154
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22 kDa
Uniprot ID	Q9H7X0
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse liver lys...