



Product Datasheet

41885 antibody (orb18375)


Description	Goat polyclonal antibody to 41885
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Septin 3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 40kDa band seen in Human Brain lysate. Recommended for use at 0.5-2 µg/ml.
Protein Sequence	SELVPEPRPKPA
Clonality	Polyclonal
MW	40.7; 40.1
Entrez	55964
Dilution Range	ELISA: 1:32000, WB: 0.5-2 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



orb18375 staining
(1ug/ml) of Human
Brai...



Immunohistochemical
staining of Human
Co...



Immunohistochemical
staining of Human
Te...