

Product Datasheet

UBE2V2 antibody (orb18381)

Description	Goat polyclonal antibody to UBE2V2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	UEV2V2 / MMS2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 20kDa band seen in human brain lysate. Recommended for use at 0.2-1 µg/ml.
Protein Sequence	AVSTGVKVP RN F
Clonality	Polyclonal
MW	16.4; 42.2
Entrez	7336
Dilution Range	ELISA: 1:32000, WB: 0.2-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human brain lys...