



Product Datasheet

DAP3 antibody (orb18383)


Description	Goat polyclonal antibody to DAP3
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Target	DAP3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 40kDa band observed in lysates of cell line HeLa (calculated MW of 45.6kDa according to NP_387506.1). In transfected HEK293 transiently expressing DAP3 a band of approx. 48kDa is observed. This band is not observed in the non-transfected HEK293. Recommended concentration: 1-3 µg/ml.
Protein Sequence	NPSLLERHCAYL
Clonality	Polyclonal
MW	45.6; 41.7; 41.0
Entrez	7818
Dilution Range	ELISA: 1:64000, WB: 0.1-1 µg/ml, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human To...



Immunohistochemical staining of HeLa lys...



Western blot analysis of HeLa lysate usi...