


Product Datasheet

TFEC antibody (orb18388)

Description	Goat polyclonal antibody to TFEC
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TFEC
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 40-45kDa band seen in Human Testis lysate. Recommended for use at 0.5-2 µg/ml.
Protein Sequence	TLDHQIINPTLK
Clonality	Polyclonal
MW	38.8; 35.5
Entrez	22797
Dilution Range	ELISA: 1:32000, WB: 0.5-2 µg/ml
Expiration Date	12 months from date of receipt.



orb18388
staining
(1ug/ml) of
Human Test...