



Product Datasheet

RGS1 antibody (orb18390)


Description	Goat polyclonal antibody to RGS1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	RGS1 / 1R20
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 25kDa band seen in HepG2 lysate. Recommended for use at 0.5-2 µg/ml.
Protein Sequence	NLLNDLQANSLK
Clonality	Polyclonal
MW	23.9
Entrez	5996
Dilution Range	ELISA: 1:32000, WB: 0.5-2 µg/ml, IHC-P: 2.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 lysate us...



Immunohistochemical staining of Human Th...



Immunohistochemical staining of Human Sm...