

Product Datasheet

KLF16 antibody (orb18392)

Description	Goat polyclonal antibody to KLF16
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	KLF16 / DRRF
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx.24kDa band observed in Human Heart lysates (calculated MW of 25.4kDa according to NP_114124.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	PAPSPAPSPAPAGL
Clonality	Polyclonal
MW	25.4
Entrez	83855
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb18392
(0.3ug/ml)
staining of
Human He...