


Product Datasheet

GRIP1 antibody (orb183922)

Description	GRIP1 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human GRIP1 (1-100/1128 aa)
Target	GRIP1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	122 kDa
Uniprot ID	Q9Y3R0
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse brain tis...