

Product Datasheet

EIF3J1 siRNA (Mouse) (orb1839336)

Description siRNA to inhibit EIF3J1 expression using RNA interference

Reactivity Mouse

Tested Applications RNA-I

Form/Appearance Lyophilized powder

Storage Shipped at 4 °C. Store at -20 °C for one year.

Note For research use only

Application notes We recommends transfection with 10 nM - 100 nM siRNA 48 to 72 hours prior to cell lysis. Before resuspending, briefly centrifuge the tube to ensure the lyophilized siRNA is at the bottom of the tube. Resuspend the siRNA oligos to an appropriate concentration with DEPC water. For example, resuspend one tube of 5 nmol siRNA oligo in 250 µl of DEPC water to get a final concentration of 20 µM.

Purity > 97%

Source Synthetic

Uniprot ID [Q3UGC7](#)

Entrez [78655](#)

Expiration Date 12 months from date of receipt.