



Product Datasheet

KLF8 antibody (orb18394)

Description	Goat polyclonal antibody to KLF8
Species/Host	Goat
Reactivity	Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	KLF8 / BKLF3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 50kDa band observed in Human Kidney lysates (calculated MW of 39.3 according to NP_009181.2). In glioblastomas this protein was reported to run at approx. 60kDa (Ding et al. J Biol Chem. 2005 Feb 25;280(8):6802-15; PMID: 15557280). Recommended for use at 1-3 µg/ml.
Protein Sequence	VDMDKLINNLEVQ
Clonality	Polyclonal
MW	39.3; 27.2
Entrez	11279
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



orb18394 staining (2ug/ml) of Human Kidn...



orb18394 (3.75µg/ml) staining of paraff...