


Product Datasheet

ARL4A antibody (orb18395)

Description	Goat polyclonal antibody to ARL4A
Species/Host	Goat
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ARL4A
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:1000.WB: Approx 23kDa band seen in Human Testis lysate. Recommended for use at 1-3 µg/ml.
Protein Sequence	GNGLSDQTSILSN
Clonality	Polyclonal
MW	22.6
Entrez	10124
Dilution Range	ELISA: 1:1000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb18395
staining
(2ug/ml) of
Human Test...