



Product Datasheet

PAX5 antibody (orb18400)


Description	Goat polyclonal antibody to PAX5
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Target	PAX5 / BSAP
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 50kDa band seen in U937 lysate (predicted MW 45kDa). Recommended for use at 2-4 µg/ml.
Protein Sequence	DLEKNYTPRTSR
Clonality	Polyclonal
MW	42.1
Entrez	5079
Dilution Range	ELISA: 1:64000, WB: 0.1-0.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of U937 lysate using PAX5 antibody.



Western blot analysis of Human Lymph lysate using PAX5 antibody.



Western blot analysis of Mouse and Rat Spleen lysates using PAX5 antibody.