



Product Datasheet

Histone H4 (tri methyl K20) antibody (orb184011)

Description	Histone H4 (tri methyl K20) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised methylpeptide derived from human Histone H4 around the methylation site of tri methyl Lys20 (HR(tri methyl K)V L)
Target	Histone H4 (tri methyl K20)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	11 kDa
Uniprot ID	P62805
Dilution Range	IHC-P=1:100-500, IHC-F=1:100-500, ICC=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
Mouse liver
using Hist...



IHC-P
analysis of
human lung
carcinoma
u...