



Product Datasheet

Hemoglobin subunit delta antibody (orb184022)

Description	Hemoglobin subunit delta antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Hemoglobin subunit delta (2-80/147 aa)
Target	Hemoglobin subunit delta
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	16 kDa
Uniprot ID	P02042
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human heart tissue.



Immunohistochemical staining of human kidney tissue.