


Product Datasheet

CABP1 antibody (orb18404)

Description	Goat polyclonal to CABP1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CABP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognise both Human isoforms.ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 20kDa band seen in human brain lysate (predicted MW 19kDa). Recommended for use at 0.5-2 µg/ml.
Protein Sequence	GNCVKYPLRNLS
Clonality	Polyclonal
MW	26.0; 19.4
Entrez	9478
Dilution Range	ELISA: 1:32000, WB: 0.5-2 µg/ml
Expiration Date	12 months from date of receipt.



orb18404
staining
(1ug/ml) of
human brai...