



Product Datasheet

POU2AF1 antibody (orb18408)

Description	Goat polyclonal antibody to POU2AF1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	BOB1 / OCAB / OBF1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 30kDa band seen in mouse spleen lysate (Predicted MW of 28kDa according to NP_006226). Recommended for use at 1-3 µg/ml.
Protein Sequence	DAYALNHTLSVEGF
Clonality	Polyclonal
MW	27.4
Entrez	5450
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse spleen ly...



Western blot analysis of HEK293 lysate u...