

Product Datasheet

EZH1 antibody (orb18414)

Description	Goat polyclonal antibody to EZH1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	EZH1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx 85kDa band observed in lysate of cell line MOLT4 (calculated MW of 85.3kDa according to NP_001982.2). Recommended concentration: 2-6 µg/ml.
Protein Sequence	EIPNPPTSKCITY
Clonality	Polyclonal
MW	85.3
Entrez	2145
Dilution Range	ELISA: 1:2000, WB: 2-6 µg/ml
Expiration Date	12 months from date of receipt.



orb18414
(2ug/ml)
staining of
lysate of ...