

Product Datasheet

HLA Class II DR4 antibody (orb184190)

Description	HLA Class II DR4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human HLA Class II DR4 (61-160/254 aa)
Target	HLA Class II DR4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	26 kDa
Uniprot ID	P01903
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.