



Product Datasheet

NOS1AP antibody (orb18425)

Description	Goat polyclonal antibody to NOS1AP
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	CAPON / NOS1AP
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.IHC: In paraffin embedded human liver shows nuclear staining in select hepatocytes. Recommended concentration, 5ug/ml.WB: No signal obtained yet but low background observed in human brain and human testis lysates at upto 4 µg/ml.
Protein Sequence	PSKTKYNLVDDGH
Clonality	Polyclonal
MW	56.1; 55.6
Entrez	9722
Dilution Range	ELISA: 1:1000, IHC-P: 4-6µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human li...



Immunohistochemical staining of Human Br...