



Product Datasheet

HTRA4 antibody (orb184261)


Description	HTRA4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human HTRA4 (251-350/476 aa)
Target	HTRA4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49 kDa
Uniprot ID	P83105
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



WB analysis of
(human/cell line)
DU145 I...



Immunohistochemical
staining of rat panc...



Western blot analysis
of (mouse/tissue) ...