



Product Datasheet

ARIH1 antibody (orb18430)


Description	Goat polyclonal antibody to ARIH1
Species/Host	Goat
Reactivity	Bovine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ARIH1 / HHARI
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 65kDa band seen in U937 lysate [Predicted MW of approx. 61kDa according to NP_005735]. Recommended for use at 0.2-1 µg/ml.
Protein Sequence	HEGYEKDLWEYIED
Clonality	Polyclonal
MW	64.1
Entrez	25820
Dilution Range	ELISA: 1:32000, WB: 0.2-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb18430
(5µg/ml)
staining of
paraffin ...



orb18430
staining
(0.5µg/ml) of
U937 ly...



Western blot
analysis of
U937 lysate
usi...