



Product Datasheet

UBE2C antibody (orb18433)

Description	Goat polyclonal antibody to UBE2C
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IP, WB
Target	UBE2C / UBCH10
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approximately 20KDa band was detected in mitotic (i.e Nocodazole arrested) HeLa whole cell lysates Recommended for use at 1-3 µg/ml. Immunoprecipitation: 20KDa band precipitated from mitotic HeLa whole cell lysates using protein-G dynabeads. Immunofluorescence: Stained mitotic and late G2 HeLa cells, with very faint or no signal in G1 and S phase cells (fixed using paraformaldehyde, primary concentration 2.5 µg/ml, using FITC anti-Goat)
Protein Sequence	QETYSKQVTSQEP
Clonality	Polyclonal
MW	19.7; 17.4; 16.7; 15.8; 14.3
Entrez	11065
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



Immunofluorescence analysis of MCF7 cell...



Western blot analysis of HEK293 (A) and ...