


Product Datasheet

PRDM16 antibody (orb18435)

Description	Goat polyclonal antibody to PRDM16
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, FC
Target	PRDM16 / MEL1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: No signal obtained yet but low background observed in Raji, Daudi, MOLT-4 and human brain lysates at up to 2µg/ml.
Protein Sequence	TSESGAFHPINHL
Clonality	Polyclonal
MW	140; 138
Entrez	63976
Dilution Range	ELISA: 1:32000, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
HEK293
cells...