


Product Datasheet

Insulin Receptor (phospho-Tyr999) antibody (orb184367)


Description

Insulin Receptor (phospho-Tyr999) antibody


Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Insulin Receptor around the phosphorylation site of Tyr999 (PE(p-Y)LS)
Target	Phospho-Insulin Receptor (Tyr999)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	affinity purified by Protein A
MW	70/152kDa
CAS Number	3
Uniprot ID	P06213
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=1µg/Test, IHC-F=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Blank
contro...



Western blot
analysis of
Lovo
(Human)Cel...



IHC-P image
of rat brain
tissue using
US...