


Product Datasheet

SNX18 antibody (orb18437)

Description	Goat polyclonal antibody to SNX18
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	SNX18 / SNAG1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: No signal obtained yet but low background observed in U937, K562 and Human Testis lysates at upto 2 µg/ml.
Protein Sequence	KVQRVPLMTVLSF
Clonality	Polyclonal
MW	72.2
Entrez	112574
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 lysate using...