



Product Datasheet

IL-12 beta antibody (orb184406)


Description	IL-12 beta antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human IL-12 beta (225-328/328 aa)
Target	IL-12 beta
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	affinity purified by Protein A
MW	35kDa
Uniprot ID	P29460
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



IHC-P image of human lung cancer tissue ...



IHC-P image of human lung carcinoma tiss...