

Product Datasheet

IL23 P19 antibody (orb184437)

Description	IL23 P19 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human IL23 P19 (51-150/189 aa)
Target	IL23 P19
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	19 kDa
Uniprot ID	Q9NPF7
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.