

Product Datasheet

IP3 receptor (phospho-Ser1598) antibody (orb184461)

Description

IP3 receptor (phospho-Ser1598)

Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human IP3 receptor around the phosphorylation site of Ser1598 (RD(p-S)V L)
Target	phospho-IP3 receptor (Ser1598)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	314 kDa
Uniprot ID	Q14643
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.