


Product Datasheet

Josephin-2 antibody (orb184482)

Description	Josephin-2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Josephin-2 (91-188/188 aa)
Target	Josephin-2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21 kDa
Uniprot ID	Q8TAC2
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Sample:
Sp2/0
(Mouse) Cell
Lysate at 30

...