


Product Datasheet

KCNA3 antibody (orb184492)

Description	KCNA3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human KCNA3 (485-575/575 aa)
Target	KCNA3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	63(Hu)/58(Mo, Ra) kDa
Uniprot ID	P22001
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 lysates u...