


Product Datasheet

KCC2 antibody (orb184496)

Description	KCC2 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human KCC2 (851-950/1139 aa)
Target	KCC2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	123 kDa
Uniprot ID	Q9H2X9
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...