



Product Datasheet

Annexin A11 antibody (orb18454)


Description	Goat polyclonal to Annexin A11
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Annexin A11
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	NP_001148.1; NP_665875.1 and NP_665876.1 are variants that encode the same protein.ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded Human Placenta shows strong staining of Basement Membrane in trophoblasts. Recommended concentration, 3-5ug/ml.WB: Approx 55kDa band seen in HeLa lysate (Predicted MW of approx. 58kDa according to NP_001148.1). Recommended for use at 0.2-1 µg/ml.
Protein Sequence	YPGYPPPPGGYPP
Clonality	Polyclonal
MW	54.4
Entrez	311
Dilution Range	ELISA: 1:32000, WB: 0.01-0.03 µg/ml, IHC-P: 3-5µg/ml
Expiration Date	12 months from date of receipt.



orb18454 staining (0.2ug/ml) of HeLa lys...



orb18454 (3.8ug/ml) staining of paraffin...



orb18454 (0.01µg/ml) staining of HeLa (...)