


Product Datasheet

SNX26 antibody (orb18457)

Description	Goat polyclonal antibody to SNX26
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	SNX26 / TCGAP (C Terminus)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Protein Sequence	SWSLHSEGGQTRSYC
Clonality	Polyclonal
MW	121
Entrez	115703
Dilution Range	ELISA: 1:8000, WB: 1 µg/ml, IHC-P: 2.5µg/ml
Expiration Date	12 months from date of receipt.



orb18457
(2.5ug/ml)
staining of
paraffin...