

Product Datasheet

Peg12 antibody (orb18465)

Description	Goat polyclonal to Peg12
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Peg12 / Frat3 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: No signal obtained yet but low background observed in Human Brain, Human Spleen, Jurkat and NIH-3T3 lysates at upto 3 µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	SLAANNPGGHSTNL
Clonality	Polyclonal
MW	29.5
Dilution Range	ELISA: 1:2000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.