


Product Datasheet

ALS2 antibody (orb18466)

Description	Goat polyclonal antibody to ALS2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	Alsin / ALS2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.IHC: In paraffin embedded Human Cortex shows pixulate cytoplasm staining in some neuronal cells. Recommended concentration, 3-5ug/ml.WB: No signal obtained yet but low background observed in Human Brain extracts at upto 1 µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	LKACYYQIQREKLN
Clonality	Polyclonal
MW	184
Entrez	57679
Dilution Range	ELISA: 1:8000, WB: 1 µg/ml, IHC-P: 3-5µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Co...