

Product Datasheet

TRIM69 antibody (orb18468)

Description

Goat polyclonal antibody to TRIM69

Species/Host

Goat

Reactivity

Bovine, Canine, Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications
Target

RNF36 / TRIM69

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:32000. WB: Preliminary experiments gave an approx 20-25kDa band in Human Testis lysates after 1 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 57.4kDa according to NP_892030.3 and of 39.3kDa according to NP_542783.2. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

GENKEPLHILHPQ

Clonality

Polyclonal

MW

57.4; 39.3

Entrez
[140691](#)
Dilution Range

ELISA: 1:32000, WB: 0.5-2 µg/ml, IHC-P: 2.5 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of Human Testis lysates



Western blot analysis of Mouse Testis lysates