


Product Datasheet

TRAK2 antibody (orb18470)

Description	Goat polyclonal antibody to TRAK2
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Target	ALS2CR3 / GRIF1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: No signals obtained yet but low background observed in Human Heart and Human Brain extracts. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any splice variants/modified forms been reported
Protein Sequence	CTSSPKMGVLKED
Clonality	Polyclonal
MW	101
Entrez	66008
Dilution Range	ELISA: 1:32000, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Flow
Cytometry
analysis of
HeLa cells
of...