

Product Datasheet

KIAA0753 antibody (orb184708)

Description	KIAA0753 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human KIAA0753 (501-600/967 aa)
Target	KIAA0753
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	109 kDa
Uniprot ID	Q2KHM9
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.