


Product Datasheet

TRIM3 antibody (orb18474)


Description	Goat polyclonal antibody to TRIM3
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	BERP / RNF22
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 85kDa band seen in Mouse Brain and Human Brain (Cerebellum) lysate [Predicted MW of approx. 80.8kDa according to NP_006449.2]. Recommended for use at 0.1-1 µg/ml.Experiment Notes: Please note that the HEK cell experiments were performed using a previous batch.
Protein Sequence	AKREDSPGPEVQP
Clonality	Polyclonal
MW	80.8
Entrez	10612
Dilution Range	ELISA: 1:64000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Brain lysate...



Western blot analysis of HEK293 lysate u...



Western blot analysis of using TRIM3 ant...