

Product Datasheet

KIAA2018 antibody (orb184757)

Description	KIAA2018 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human KIAA2018 (301-400/2245 aa)
Target	KIAA2018
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	241 kDa
Uniprot ID	Q68DE3
Expiration Date	12 months from date of receipt.