

Product Datasheet

CECR2 antibody (orb18494)

Description

Goat polyclonal to CECR2

Species/Host

Goat

Reactivity

Canine, Human

Conjugation

Unconjugated

Tested

ELISA, FC, IF

Applications

Target

CECR2

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

The peptide sequence used (from AAK15343.1) differs by one residue from the RefSeq entry (XP_931107.2). It is not known whether this difference is due to a sequencing error or polymorphism but the antibody is expected to recognize both sequences. ELISA: Peptide ELISA: antibody detection limit dilution 1:32000. WB: No signal obtained yet but low background observed in human testis and 293 lysates at up to 1 µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates

Protein Sequence

PPHKPPTLPLDQS

Clonality

Polyclonal

MW

164;162;143

Entrez


27443

Dilution Range


ELISA: 1:32000, IF/ICC: 10µg/ml, FACS: 10ug/ml

Expiration Date

12 months from date of receipt.



Flow cytometric analysis of paraformalde...



Immunofluorescence analysis of paraforma...